User Manual

English

TERMAKS

TERB9025

TERB9051

TERB9130

TERB9420

The original document B9000
of the manufacturer is obligatory.

User's Manual

Incubators


Models

B 9025 - B 9051 - B 9130 - B 9420

Content /Innhold	Page/Side
Basic Functions / <i>Grunnfunksjoner</i>	2 / 10
Installation / <i>Innstallasjon</i>	3 / 11
Main Switch and External Connections / <i>Hovedbryter og tilkoblinger</i>	3 / 11
Display Only Frames, Overview / <i>Displayoversikt, aktuelle verdier</i>	4 / 12
Data Logging / <i>Datalogging</i>	4 / 12
Settings Menu / <i>Innstillingsmeny</i>	5 / 13
Timers / <i>Timere</i>	6 / 14
Sterilizing Cycle / <i>Steriliseringssyklus</i>	6 / 14
The Fan / <i>Viften</i>	7 / 15
The Air Valve / <i>Spjeldet</i>	7 / 15
Standby / <i>Standby</i>	7 / 15
Pin Code / <i>Adgangskode</i>	8 / 16
Alarms / <i>Alarmer</i>	8 / 16
Safety Thermostat / <i>Sikkerhetstermostat</i>	9 / 17
Calibration / <i>Kalibrering</i>	9 / 17

Basic Functions

This page gives you all information needed to operate the chamber. Details are given on the following pages.

How to Change a Setting

- Go to Display **Settings Mode** by pressing the **Set Key**. Set LED turns on and explaining text appears in left display
- Find desired parameter in the **Settings Menu** with **Up/Down keys**
- Go to **Change Mode** by pressing **Set Key** for 3 sec. LED turns flashing and explaining text appears in left display
- **Adjust** the setting with **Up/Down Keys**
- **Save** the new setting by pressing the **Enter Key**

How to Operate the Standard Timer

- Enter **CycleTime** (hours and minutes) as described above
- **Start** timer by pressing **Timer Key** until the LED turns on
- **Heating Up** is displayed until the Set temperature has been reached
- **Time Left** (Hours Minutes) will be displayed when the set temperature has been reached.


When Cycle is Complete:


- **Cycle Complete** appears in the display
- The LED starts flashing
- Heating turns off
- Air valve closes (If assembled, normally not)
- Fan stops (8 minutes delayed)
- The cycle can be **repeated** by pressing **Timer Key** until LED turns on

Installation

- Unpack and place shelves in desired position
- Place the chamber with at least 10 cm free air on all sides
- Do not connect the air outlet (if assembled) to a ventilation system
- Connect the chamber to a grounded power source
- Turn it on with the **Main Switch** on the right side.
- The display turns on and shows **Wait...** for a few seconds
- **Set** the desired **Temperature, Fan Speed**, and eventually, the **Time**
- If Norwegian is preferred in the display, go to System Setup and change language
- If **Real Time Program** is preferred, choose it in the **System Setup menu**

Main Switch and External Connections


Temperature Logging

- Starts when the temperature has been reached, about 30 minutes delayed
- Used to evaluate the temperature control
- Updates every minute
- Door openings and Air Valve adjustment will affect the Max and Min values
- All zeros means that the logging has not yet started
- Cleared when Power Up, wakeup from Standby, entering new Temperature Set and when turning on the Timer

Settings Menu


Timers

In the **System Setup**, one of **two timer systems** can be chosen:

Standard Timer (default):

- Can be set to maximum 99 hours and 59 minutes
- **Starts counting** down when the **Set Temperature** is reached
- In the case of power failure, lasting more than 5 minutes, the counting will restart from the beginning
- Operated as described in the chapter **How to Operate the Standard Timer**
- The Fan and Air valve (if assembled) are controlled as illustrated below
- If the chamber shall run without any timing, the time setting should be set to 00h 00m. This will avoid accidental start of the timer



Sterilization Cycle

- A **fixed sterilization cycle** for the chamber is available , 180°C in 40 minutes
- **Standard timer** must have been chosen in the **Setup Menu**
- Behave as Standard Timer, but temperature and time can't be adjusted
- The cycle can be **started** by pressing **Up, Down and Timer keys** for 3 seconds
- The message **Sterilization** appears in the left display and the Timer LED turns on
- **Heating Up** is displayed until the 180°C has been reached
- **Time Left** (Minutes) will be displayed when the set temperature has been reached.
- When finished, **Cycle Complete** will be displayed and the LED starts flashing
- Return to normal operation, timer or continuous, by pressing Timer key for 3 sec

Real Time Program:

- Optional, can be chosen in the System Setup Menu
- Makes it possible to operate the chamber at two different temperatures
- The extra temperature is named **Temperature B** in the menus
- Setting Temperature B to 00.0°C means that heater and fan are turned off in that period. (Sleep Mode)
- **Start- and end Time** for Temperature B can be set to any time during the **week or on daily basis**
- When choosing **xx** instead of a specified day, the change will occur daily at the given time
- The **Real Time Clock** is backed up for about 10 minutes if power fails. A warning is given if adjustment is needed.

The Fan

- Adjustable speed, 1 to 10
- When heating up, settings lower than 6 will be ignored
- Stops for about 30 seconds when the door is opened while running at set temperature. This function is disabled for 4 minutes after last door opening or after last Air Valve adjustment
- May reduce speed, or even stop, if operating at a temperature close to the ambient temperature

The Air Valve (if assembled)

- Adjustable opening in steps of 10% (00 = closed, 100 = open)
- Movement may be 10 sec delayed
- Should be open when used for drying purposes
- Closes automatically when heating up, when a timer cycle is completed or in Sleep mode (Real Time Program)

Standby

- Preferred way to shut down the chamber. Keeps the real time clock running
- Turns off Heaters, Fan and Display
- Press the Standby Key for about 3 seconds until the backlight in the display turns off
- The Standby LED turns on when the key is released.
- The chamber can be restarted by pressing the Standby Key for 3 seconds
- The right display will show Wait... for a moment during initialization

Pin Code

- Optional system to stop unwanted settings from unauthorized persons
- A **personal pin code** can be entered during the first 30 seconds after turning the chamber on with the main switch or waking up from Standby Mode. Actual code is displayed after pressing the Set key in that period. A unique code can be chosen in the range from 001 to 999.
- 000 means deactivated Pin Code System.
- When activated, **Pin Code?** will be displayed in the left display.
- **Actual code must be set by Up/Down Keys**, followed by pressing the **Enter Key**
- The system will remain opened for 30 seconds after the last key operation.
- Alarms can be acknowledged by the Standby Key even if the keys are locked
- If assembled with door lock system, the door will be locked as the keys are locked.

Alarms

- A number of alarm situations may occur, listed below
- Actual warning is given in the lower line of the left display
- The backlight starts flashing and the buzzer starts sounding whenever an alarm occur
- The flashing and the sound can be stopped by pressing the Standby Key

Displayed Warning Left Display	Explanation	
High Temp Alarm	The temperature has drifted above the entered High Temperature Alarm Limit Turns off 0,5°C below the Alarm Limit	<p style="text-align: center;">Disabled during: Change Up/ Down, Cycle Complete and Sleep Mode</p> <p style="text-align: center;">Disabled for 4 minutes after:</p> <ul style="list-style-type: none"> • Door Opening • Valve Adjustment • Heating Up • Cooling Down
Low Temp Alarm	The temperature has dropped below the entered High Temperature Alarm Limit Turns off 0,5°C above the Alarm Limit	
High Temp Safety	The Temperature Safety Limit has been exceeded and the power to the heater is shut off	
Sensor Failure	One of the two sensors, for display or safety, is probably broken. The chamber operates with limited safety	
Adjust Real Time!	The Real Time Clock needs adjustment. Probably due to a significant power failure	
Data Error	Communication error between the panel and the controller. The chamber operates with the last valid settings, but the displayed data may be wrong	

Alarm Output

- A potential free contact is available for external use (see drawing Main Switch and External Connections)
- The output switches when an alarm situation occurs, one minute delayed
- Switches immediately at power failure
- Immediate return to normal when the alarm situation is over

Safety Thermostat

- Automatic system, preventing risky temperatures in the case of a controller failure
- Completely independent from the ordinary temperature control system
- Automatically adjusted whenever a new Temperature Set is entered
- Adjustable limit (offset from target)
- Cuts power to the heaters when the limit is exceeded
- Generates an alarm warning
- Turns to normal when the temperature has dropped 0,5°C below the limit

Calibration

Two calibration constants can be adjusted in order to bring the temperature measurement equal to the true value according to a control device:

Amplifying Factor (A)

Should be calculated if the chamber is used in a wide temperature range

Zero Point Constant (Z)

Should be calculated if display and a control device is not equal at normal working temperature

The calibration should be performed as follows:

Place the control device into the middle of the cabinet and let the cabinet stabilize at a low temperature
Note the displayed temperature, T_{dl} and the true temperature measured in the control device T_{cl} .

Example:

$$T_{dl} = 30,0^{\circ}\text{C} - T_{cl} = 30,8^{\circ}\text{C}$$

Let the cabinet stabilize at a high temperature
Note the displayed temperature, T_{dh} and the true temperature measured in the control device T_{ch}

$$T_{dh} = 80,0^{\circ}\text{C} - T_{ch} = 81,5^{\circ}\text{C}$$

Find current amplifying factor in the Settings Menu

$$A = 750$$

Calculate the new factor:

$$A_{\text{new}} = \{(T_{ch} - T_{cl}) / (T_{dh} - T_{dl})\} * A$$

Enter A_{new}

$$A_{\text{new}} = \{(81,5 - 30,8) / (80,0 - 30,0)\} * 750 = 761$$

Let the cabinet stabilize at a temperature in the middle of the normal working range

Note the displayed temperature T_d and the true temperature measured in the control device T_c

$$T_d = 40,0^{\circ}\text{C} \quad T_c = 40,8^{\circ}\text{C}$$

Find current zero point constant in the Settings Menu

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